

MiniDygestorium-350/Ex – individual stand for work with dusts and gases



II 2 G c Ex e II T3



- pre-filter,
- high-efficiency HEPA filter – class H13,
- gas absorber – a cassette with granular activated carbon,
- Ex fan placed in the lower part of the device, at the side of clean air,
- pressure control – indicating the excessive resistances of the high-efficiency filter,
- control unit (to be installed within the room, beyond the Ex hazard area).

Operational Use

The construction is an independent mobile workplace. After switching it on, the operator places the emission source on the desktop (inside the cabinet), whereby the tasks are executed in the vacuum area, that eliminates the pollution being emerged outside.

The dust pollutants are captured by the pre-filter and the high-efficiency HEPA filter. Whereas, the active carbon layer absorbs the majority of noxious chemical compounds, such as: styrene, toluene, alcohols, phenol and many others. At the point when the HEPA filter reaches the limit pollution degree, a light signal indicates the need of filter replacement.

Air is supplied into the extraction cabinet through the perforated upper wall and the holes for hands (in the front). The polluted air is expelled through the perforated outlet, located underneath the device.

Maintenance consists in:

- periodical replacement of the HEPA filter – as signalled by the lamp,
- periodical replacement of the cassette with active carbon – depending on organoleptic evaluation of operator,
- periodical replacement of the pre-filter.

CAUTION:

Absorption efficiency of the active carbon for various vapours and gases is listed on the next page.

Purpose

MiniDygestorium-350/Ex has been developed for purifying the air of the gaseous contaminations, emitted in small amounts, in chemical laboratories, biological-, analytical-, scientific facilities, research labs, health service units, in chemical ateliers in schools and in numerous other places, where noxious gases and vapours arise, which endanger our health.

MiniDygestorium-350/Ex eliminates the expansion possibility of the pollutants within the room. The appliance can be used in areas of explosion hazard, where explosive atmosphere is likely to occur.

Structure

The device consists of following elements:



- cabinet fume hood – a glass extraction chamber made of acid-proof steel, with two holes for operator's hands, due to which various operations can be carried out on the desktop,
- housing of steel sheets – 3 segments assembled together with clasp locks,

Technical Data

Type	Part No.	Maximum volume flow [m³/h]	Maximum vacuum [Pa]	Motor rate [W]	Supply voltage [V/Hz]	Acoustic pressure level [dB(A)]*	Weight [kg]
MiniDygestorium-350/Ex	888D01	350	220	120	3x400	48	98


* Noise level has been measured at a distance of 1 metre (from the device).

Replaceable Parts


High-efficiency HEPA filter							Cassette with activated carbon					
	Type	Part No.	Weight [kg]	Dimensions AxB xH [mm]	Class	Filtration material		Type	Part No.	Weight [kg]	Dimensions AxB xH [mm]	Remarks
	FW-MD-350/Ex	838W03	15	535x535 x292	H13	Hydrophobic glass paper 99,95%		WA-ECO-20	838K98	24*	534x534 x155	The cassette is made of cardboard and plywood

*Weight of the activated carbon ~20 kg

Pre-filter

	Type	Part No.	Weight [kg]	Dimensions AxB xH [mm]	Class	Filtration material
	PS-MD-350	852F03	0,5	535x535 x50	G3	Glass unwoven with progressively growing density

Dimensions



Values of activated carbon absorption efficiency for various types of vapors and gases

High efficiency

ethyl acrylate - $C_5H_8O_2$
 methyl acrylate - $C_4H_6O_2$
 acrylo_nitrile - C_3H_3N
 valeraldehyde - $C_5H_{10}O$
 amyl alcohol - $C_5H_{12}O$
 butyl alcohol - $C_4H_{10}O$
 propyl alcohol - C_3H_7OH
 a_nili_ne - $C_6H_5NH_2$
 „aphta (petroleum)
 „aphta (coal tar)
 bromi_ne - Br_2
 butyl cellosolve - $C_6H_{14}O_2$
 - cellosolve - $C_4H_{10}O_2$
 - cellosolve acetate - $C_6H_{12}O_3$
 butyl chloride - C_4H_9Cl
 propyl chloride - C_3H_7Cl
 mo_nochlorobe_nze_ne - C_6H_5Cl
 chlorobe_nze_ne - C_6H_5Cl
 ethyle_ne chlorhydri_ne - C_2H_5ClO
 chloroform - $CHCl_3$
 chloro_nitropropa_ne - $C_3H_6ClNO_2$
 chloropici_ne - CCl_3NO_2
 chlorobutadiene - C_4H_5Cl
 cyclohexa_nol - $C_6H_{12}O$
 cyclohexa_no_ne - $C_6H_{10}O$
 tetrachloroetha_ne - $C_2H_2Cl_4$
 tetrachloroethyle_ne - C_2Cl_4
 carbo_n tetrachloride - CCl_4
 deca_ne - $C_{10}H_{22}$
 dioxa_ne - $C_4H_8O_2$
 dibromometha_ne - CH_2Br_2
 ethyle_ne dichloride - $C_2H_4Cl_2$
 dichlorobe_nze_ne - $C_6H_4Cl_2$
 dichloroetha_ne - $C_2H_4Cl_2$
 dichloroethyle_ne - $C_2H_2Cl_2$
 dichloro_nitroetha_ne - $CH_3CCl_2NO_2$
 dichloropropa_ne - $C_3H_5Cl_2$
 dimethyla_nili_ne - $C_6H_{11}N$
 amyl ether - $C_{10}H_{22}O$
 butyl ether - $C_8H_{18}O$
 dichloroethyl ether - $C_4H_8Cl_2O$
 isopropyl ether - $C_6H_{14}O$
 propyl ether - $C_6H_{14}O$
 ethyl be_nze_ne - C_8H_{10}
 phe_nol - C_6H_5O
 hepta_ne - C_7H_{16}
 heptyle_ne - C_7H_{14}
 i_ndole - C_8H_7N
 isophoro_ne - $C_9H_{14}O$
 iodi_ne - I
 iodoform - CHI_3
 camphor - $C_{10}H_{16}O$
 diethyl keto_ne - $C_5H_{10}O$

dipropyl keto_ne - $C_7H_{14}O$
 methyl butyl keto_ne - $C_6H_{12}O$
 methyl isobutyl keto_ne - $C_6H_{12}O$
 methyl ethyl keto_ne - C_4H_8O
 creosole - $C_8H_{10}O_2$
 cresol - C_7H_8O
 croto_naldehyde - C_4H_6O
 ethyl silicate - $C_8H_{20}O_4Si$
 acrylic acid - $C_3H_4O_2$
 caprylic acid - $C_8H_{16}O_2$
 butyric acid - $C_4H_8O_2$
 lactic acid - $C_3H_6O_3$
 uric acid - $C_5H_4N_4O_3$
 acetic acid - CH_3COOH
 propio_nic acid - $C_3H_6O_2$
 valeric acid - $C_5H_{10}O_2$
 me_nthol - $C_{10}H_{20}O$
 ethyl mercapta_ne - C_2H_5S
 propyl mercapta_ne - C_3H_7S
 - methyl cellosolve - $C_3H_8O_2$
 - methyl cellosolve acetate - $C_5H_{10}O_3$
 methylcyclohexa_ne - C_7H_{14}
 methylcyclohexa_nol - $C_7H_{14}O$
 urea - CH_4N_2O
 kero_nse_ne
 „icoty_ne - $C_{10}H_{14}N_2$
 „itrobe_nze_ne - $C_6H_5NO_2$
 „itroetha_ne - $C_2H_5NO_2$
 „itroglyceri_ne - $C_3H_5N_3O_9$
 „itropropa_ne - $C_3H_7NO_2$
 „itrotolue_ne - $C_7H_7NO_2$
 „o_na_ne - C_9H_{20}
 amyl acetate - $C_7H_{14}O_2$
 butyl acetate - $C_6H_{12}O_2$
 ethyl acetate - $C_4H_8O_2$
 isopropyl acetate - $C_5H_{10}O_2$
 propyl acetate - $C_5H_{10}O_2$
 octale_ne - $C_{12}H_{26}Cl_6$
 octa_ne - C_8H_{18}
 putresci_ne - $C_4H_{12}N_2$
 ozo_ne - O_3
 paradichlorobe_nze_ne - $C_6H_4Cl_2$
 - pe_nta_no_ne - $C_5H_{10}O$
 perchloroethyle_ne - C_2Cl_4
 pyridi_ne - C_5H_5N
 dimethylsulphate - $C_2H_6O_4S$
 skatole - C_9H_7N
 styre_ne mo_nomer - C_8H_8
 turpe_nti_ne - $C_{10}H_{16}$
 mesityl oxide - $C_6H_{10}O$
 tolu_ne - C_7H_8
 toluidi_ne - C_7H_9N
 trichloroethyle_ne - C_2HCl_3

Average efficiency

aceto_ne - C_3H_6O
 acetyle_ne - C_3H_2
 acrolei_ne - C_3H_4O
 butyraldehyde - C_4H_8O
 ethyl alcohol - C_2H_5OH
 methyl alcohol - CH_3OH
 be_nze_ne - C_6H_6
 ethyl bromide - C_2H_5Br
 methyl bromide - CH_3Br
 butadie_ne - C_4H_6
 chlori_ne - Cl_2
 ethyl chloride - C_2H_5Cl
 vi_nyl chloride - C_2H_3Cl
 cyclohexe_ne - C_6H_{10}
 dichlorodifluorometha_ne - CCl_2F_2
 diethyl ami_ne - $C_4H_{11}N$
 carbo_n disulphide - CS_2
 ether - $C_4H_{10}O$
 ethyl ether - $C_4H_{10}O$
 ethyl ami_ne - C_2H_7N
 fluorotrichlorometha_ne - CCl_3F
 phosge_ne - $COCl_2$
 a_naesthetics
 hexa_ne - C_6H_{14}
 hexyle_ne - C_6H_{12}
 hexy_ne - C_6H_{10}
 isopre_ne - C_5H_8
 hydroge_n iodide - HI
 xyle_ne - C_8H_{10}
 formic acid - $HCOOH$
 methyl mercapta_ne - CH_3SH
 ethyl formate - $C_3H_6O_2$
 methyl formate - $C_2H_4O_2$
 „itrometha_ne - CH_3NO_2
 methyl acetate - $C_3H_6O_2$
 pe_nta_ne - C_5H_{12}
 pe_ntyle_ne - C_5H_8
 pe_nty_ne - C_5H_6
 propio_na_ndehyde - C_3H_6O
 ethyle_ne oxide - C_2H_4O
 carbo_n mo_noxide - CO

Low efficiency

acetaldehyde - C_2H_4O
 ammo_nia - NH_3
 hydroge_n bromide - HBr
 buta_ne - C_4H_{10}
 buta_no_ne - C_4H_8O
 butyle_ne - C_4H_8
 buty_ne - C_4H_6
 methyl chloride - CH_3Cl
 hydroge_n chloride - HCl
 hydroge_n cya_nide - HCN
 „itroge_n dioxide - NO_2
 sulphur dioxide - SO_2
 hydroge_n fluoride - HF
 formaldehyde - CH_2O
 propa_ne - C_3H_8
 propyle_ne - C_3H_6
 propy_ne - C_3H_4
 hydroge_n sele_nide - H_2Se
 hydroge_n sulphide - H_2S
 sulphur trioxide - SO_3