

MiniDygestorium-350 – individual stand for work with dusts and gases



- carbon,
- radial fan with enclosure,
- pressure control that signalizes excessive resistance on the high efficiency filter,
- control unit.

Operational Use

The device is an independent, mobile workstation. After switching it on, the operator places the source of emission on the desktop inside the cabinet hood and performs the required tasks inside the negative pressure zone, which eliminates the possibility of polluting the room with gases.

Dust contaminants are then filtered by the high efficiency HEPA filter, while the case with activated carbon absorbs most of the harmful chemical compounds, such as styrene monomer, toluene, alcohols, phenol and other. When the threshold value of HEPA filter clogging is reached, a light signal will be activated to inform the operator about the need to replace the filter.

Air gets inside the cabinet hood through the perforated upper wall of the extraction chamber and through the openings for hands located in the front wall. It is removed through the perforated outlet located beneath the device.

Operation of the device requires:

- periodic replacement of the HEPA filter – the need to replace the filter is signalized by a red lamp,
- periodic replacement of the case with activated carbon – the need to replace the case is determined on the basis of a visual inspection of the case,
- periodic replacement of the Paint-Stop filter.

NOTE:

Absorption efficiency of activated carbon differs depending on the type of vapor or gas. Referential values are given on the following page.

Purpose

MiniDygestorium-350 was designed to purify air and absorb any waste gases produced in small amounts in chemical, biological or analytical laboratories, research or healthcare facilities, school labs and other places that produce where harmful gases or vapors are produced.

MiniDygestorium-350 makes it impossible for waste gases to spread across the room. The device must not be used in hazardous areas with explosive atmospheres.

Structure

The device is composed of the following elements:

- cabinet hood – glass-covered exhaustion chamber made of acid-proof steel with two openings for hands allowing for the performance of various tasks inside the device,
- enclosure made of steel sheets – 3 segments secured with clasps,
- Paint-Stop filter,
- High efficiency HEPA filter of H13 class,
- gas absorber in the form of a case with granular activated

Technical data

| Type | Part no. | Maximal volume flow [m ³ /h] | Maximal vacuum [Pa] | Motor rate [W] | Supply voltage [V/Hz] | Acoustic pressure level [dB(A)]* | Weight [kg] |
|---------------------|----------|---|---------------------|----------------|-----------------------|----------------------------------|-------------|
| MiniDygestorium-350 | 801O20 | 350 | 220 | 124 | 230/50 | 53 | 80 |

* Measurements were performed 1 m away from the device.

Spare parts

High efficiency HEPA filter

| | Type | Part no. | Weight [kg] | Dimensions AxB xH [mm] | Class | Filter material |
|--|-----------|----------|-------------|------------------------|-------|--|
| | FW-MD-350 | 838F98 | 3,2 | 535x535 x78 | H13 | Hydrophobic glass fiber filter paper 99.95%. |

Case with activated carbon

| | Type | Part no. | Weight [kg] | Dimensions AxB xH [mm] | Remarks |
|--|-----------|----------|-------------|------------------------|--|
| | WA-ECO-20 | 838K98 | 24* | 534x534 x155 | Case is made of cardboard and plywood. |


*Weight of activated carbon 20 kg.

Pre-filter

| | Type | Part no. | Weight [kg] | Dimensions AxB xH [mm] | Class | Filter material |
|--|-----------|----------|-------------|------------------------|-------|--|
| | PS-MD-350 | 852F03 | 0,5 | 535x535 x50 | G3 | Non woven glass fabric of progressive density. |

MiniDygestorium-350

Dimensions



Values of activated carbon absorption efficiency for various types of vapors and gases

High efficiency

ethyl acrylate - C₅H₈O₂
 methyl acrylate - C₆H₈O₂
 acrylo_nitrile - C₃H_{3.5}N
 valeraldehyde - C₅H₁₀O
 amyl alcohol - C₅H₁₂O
 butyl alcohol - C₄H₁₀O
 propyl alcohol - C₃H₇OH
 a_nili_ne - C₈H₅NH₂
 a_naphta (petroleum)
 a_naphta (coal tar)
 bromi_ne - Br₂
 butyl cellosolve - C₆H₁₄O₂
 - cellosolve - C₄H₁₀O₂
 - cellosolve acetate - C₆H₁₂O₃
 butyl chloride - C₄H₉Cl
 propyl chloride - C₃H₇Cl
 mo_nochlorobe_nze_ne - C₆H₅Cl
 chlorobe_nze_ne - C₆H₅Cl
 ethyle_ne chlorhydri_ne - C₂H₅ClO
 chloroform - CHCl₃
 chloro_nitropropa_ne - C₃H₆ClNO₂
 chloropicri_ne - CCl₃NO₂
 chlorobutadie_ne - C₄H₅Cl
 cyclohexa_nol - C₆H₁₂O
 cyclohexa_no_ne - C₆H₁₀O
 tetrachloroetha_ne - C₂H₂Cl₄
 tetrachloroethyle_ne - C₂Cl₄
 carbo_ne tetrachloride - CCl₄
 deca_ne - C₁₀H₂₂
 dioxa_ne - C₄H₈O₂
 dibromometha_ne - CH₂Br₂
 ethyle_ne dichloride - C₂H₄Cl₂
 dichlorobe_nze_ne - C₆H₄Cl₂
 dichloroetha_ne - C₂H₄Cl₂
 dichloroethyle_ne - C₂H₂Cl₂
 dichloro_nitroetha_ne - CH₃CCl₂NO₂
 dichloropropa_ne - C₃H₄Cl₂
 dimethyla_nili_ne - C₈H₁₁N
 amyl ether - C₁₀H₂₂O
 butyl ether - C₉H₁₈O
 dichloroethyl ether - C₄H₈Cl₂O
 isopropyl ether - C₆H₁₄O
 propyl ether - C₆H₁₄O
 ethyl be_nze_ne - C₈H₁₀
 phe_nol - C₆H₆O
 hepta_ne - C₇H₁₆
 heptyle_ne - C₇H₁₄
 i_ndole - C₈H₇N
 isophoro_ne - C₉H₁₄O
 iodi_ne - I
 iodoform - CHI₃
 camphor - C₁₀H₁₆O
 diethyl keto_ne - C₅H₁₀O

dipropyl keto_ne - C₇H₁₄O
 methyl butyl keto_ne - C₆H₁₂O
 methyl isobutyl keto_ne - C₆H₁₂O
 methyl ethyl keto_ne - C₄H₈O
 creosole - C₈H₁₀O₂
 cresol - C₇H₈O
 croto_naldehyde - C₄H₆O
 ethyl silicate - C₈H₂₀O₄Si
 acrylic acid - C₃H₄O₂
 caprylic acid - C₈H₁₆O₂
 butyric acid - C₄H₈O₂
 lactic acid - C₃H₆O₃
 uric acid - C₅H₄N₄O₃
 acetic acid - CH₃COOH
 propio_nic acid - C₃H₆O₂
 valeric acid - C₅H₁₀O₂
 me_nthol - C₁₀H₂₀O
 ethyl mercapta_ne - C₂H₆S
 propyl mercapta_ne - C₃H₈S
 - methyl cellosolve - C₃H₈O₂
 - methyl cellosolve acetate - C₅H₁₀O₃
 methylcyclohexa_ne - C₇H₁₄
 methylcyclohexa_nol - C₇H₁₄O
 urea - CH₄N₂O
 kerose_ne
 i_nco_nty_ne - C₁₀H₁₄N₂
 i_ntrobe_nze_ne - C₆H₅NO₂
 i_ntroetha_ne - C₂H₅NO₂
 i_ntrogliceri_ne - C₃H₅N₂O₉
 i_ntropropa_ne - C₃H₇NO₂
 i_ntrotole_ne - C₇H₇NO₂
 n_on_oa_ne - C₃H₂₀
 amyl acetate - C₇H₁₄O₂
 butyl acetate - C₆H₁₂O₂
 ethyl acetate - C₄H₈O₂
 isopropyl acetate - C₅H₁₀O₂
 propyl acetate - C₅H₁₀O₂
 octale_ne - C₁₂H₈Cl₆
 octa_ne - C₈H₁₈
 putresci_ne - C₄H₁₂N₂
 ozo_ne - O₃
 paradichlorobe_nze_ne - C₆H₄Cl₂
 - pe_nta_no_ne - C₅H₁₀O
 perchloroethyle_ne - C₂Cl₄
 pyridi_ne - C₅H₅N
 dimethylsulphate - C₂H₆O₄S
 skatole - C₈H₉N
 styre_ne mo_nomer - C₈H₈
 turpe_nti_ne - C₁₀H₁₆
 mesityl oxide - C₈H₁₀O
 tole_ne - C₇H₈
 toluidi_ne - C₇H₉N
 trichloroethyle_ne - C₂HCl₃

Average efficiency

aceto_ne - C₃H₆O
 acetyle_ne - C₂H₂
 acrolei_ne - C₃H₄O
 butyraldehyde - C₄H₈O
 ethyl alcohol - C₂H₅OH
 methyl alcohol - CH₃OH
 be_nze_ne - C₆H₆
 ethyl bromide - C₂H₅Br
 methyl bromide - CH₃Br
 butadie_ne - C₄H₆
 chlori_ne - Cl₂
 ethyl chloride - C₂H₅Cl
 vi_nyl chloride - CH₂CHCl
 cyclohexe_ne - C₆H₁₀
 dichlorodifluorometha_ne - CCl₂F₂
 diethyl ami_ne - C₄H₁₁N
 carbo_ne disulphide - CS₂
 ether - C₂H₆O
 ethyl ether - C₄H₁₀O
 ethyl ami_ne - C₂H₇N
 fluorotrichlorometha_ne - CCl₃F
 phosge_ne - COCl₂
 a_naesthetics
 hexa_ne - C₆H₁₄
 hexyle_ne - C₆H₁₂
 hexy_ne - C₆H₁₀
 isopre_ne - C₅H₈
 hydroge_ne iodide - HI
 xyle_ne - C₈H₁₀
 formic acid - HCOOH
 methyl mercapta_ne - CH₃SH
 ethyl formate - C₃H₆O₂
 methyl formate - C₂H₄O₂
 i_ntrometha_ne - CH₃NO₂
 methyl acetate - C₃H₆O₂
 pe_nta_ne - C₅H₁₂
 pe_ntyle_ne - C₅H₈
 pe_nty_ne - C₅H₈
 propio_na_ndehyde - C₃H₆O
 ethyle_ne oxide - C₂H₄O
 carbo_ne mo_noxide - CO

Low efficiency

acetaldehyde - C₂H₄O
 ammo_nia - NH₃
 hydroge_ne bromide - HBr
 buta_ne - C₄H₁₀
 buta_no_ne - C₄H₈O
 butyle_ne - C₄H₈
 buty_ne - C₄H₆
 methyl chloride - CH₃Cl
 hydroge_ne chloride - HCl
 hydroge_ne cyanide - HCN
 i_ntroge_ne dioxide - NO₂
 sulphur dioxide - SO₂
 hydroge_ne fluoride - HF
 formaldehyde - CH₂O
 propa_ne - C₃H₈
 propyle_ne - C₃H₆
 propy_ne - C₃H₄
 hydroge_ne sele_nide - H₂Se
 sulphur sulphide - H₂S
 sulphur trioxide - SO₃